




UniGate i (UGi)

The UniGate i (UGi) is an advanced DICOM gateway that captures and converts non-DICOM images into the DICOM format. It connects non-DICOM compliant modalities to standard DICOM solutions and systems. UGi is available as a mobile application specifically designed for the iPad, leveraging its high-quality display and portability to enhance usability and mobility for medical professionals, in comparison to conventional network workstations.

Features

- » **Wireless and Lightweight Portability**
This lightweight device saves space for essential medical materials and equipment.
- » **Comprehensive Information Recording**
Provides border control and personnel records to manage infection risks.
- » **High-Definition Image Capture**
Supports 2K and 4K options, ensuring Full HD digital image capture quality.
- » **Precise Data Integration and Management**
Integrates seamlessly with HIS/Worklist systems and automates numbering for quick patient information construction.
- » **Intuitive Interface and User-Friendly Design**
Easy to navigate with one-click access to log records.
- » **Multi-Format File Export**
Supports image capture and dynamic video recording in formats like DICOM and JPG.

System Workflow



*iPad® is registered trademarks of Apple Inc.